

Bright & Early Postural Care:

Considerations for adaptive equipment support in high-risk infants

Series 2: Starting the supported standing journey

Introductions

Kayleah Cooper

National Clinical Training Manager
Physiotherapy BSc



Kieran Murphy

Clinical Partnerships Manager
Physiotherapy BSc, MSc



Overview

- What is Early Intervention?
- Learnings from Cerebral Palsy interventions
- Overview of standing frame options
- Indications/precautions/contraindications for standing
- Understanding the importance of carer and MDT collaboration and child centred communication
- Utilising these principles for individual, holistic, child-centred care

What is 'Early Intervention'?

- Depends on the context:
- Detection, treatment, management? Single therapy or multi – disciplinary?
- Early, or On time?



For the purpose of this webinar:

Early intervention: A programme to support postural and motor development in higher risk infants

Development: Skeletal









Image from National Library of Medicine 2023, reproduced U.S. Department of Labor, Children's Bureau. Infant Growth and Development, Ca. 1943

https://www.nlm.nih.gov/exhibition/dream-anatomy/index.html?assetPath=section3_slide23



Development: Gross Motor



Pull to standing
Cruise

8
10

Windows of achievement for six gross motor milestones



Reference: WHO Multicentre Growth Reference Study Group. WHO Motor Development Study: Windows of achievement for six gross motor development milestones. Acta Paediatrica Supplement 2006;450:86-95.

Retrieved from: https://cdn.who.int/media/docs/default-source/child-growth/child-growth-standards/indicators/motor-development-milestones/who-motor-development-study-windows-of-achievement-for-six-gross-motor-development-milestones.pdf?sfvrsn=3425c1dc_0


Play

- Play is an integral part of a child's life
- Helps build brain architecture
- Helps develop postural control and gross motor skills
- Encourages family and social interaction



Copyright **Pathways.org** retrieved from <https://pathways.org/kids-learn-play-6-stages-play-development/>

Development: Neuroplasticity



Center on the Developing Child HARVARD UNIVERSITY

SOURCE: LEVITT (2009)

www.developingchild.harvard.edu



'Infants and children have remarkable brain plasticity and aptitude for learning; taking advantage of this window of opportunity with evidence-based intervention for those with early brain injury maximizes functional outcomes and minimizes complications.'


Morgan et al. 2021

Developments in Early Diagnosis of Cerebral Palsy

Detection and Diagnosis

- Traditionally, CP has been diagnosed between the ages of 12 and 24 months.
- CP can now potentially be diagnosed as early as 3 months corrected gestational age (CGA)

“Early diagnosis begins with a medical history and involves using neuroimaging, standardized neurological, and standardized motor assessments that indicate congruent abnormal findings indicative of cerebral palsy. Clinicians should understand the importance of prompt referral to diagnostic-specific early intervention to optimize infant motor and cognitive plasticity, prevent secondary complications, and enhance caregiver well-being” (Novak et al, 2017)



(Novak et al, 2017 – accessed from American Academy for Cerebral Palsy and Developmental Medicine: <https://www.aacpdm.org/publications/care-pathways/early-detection-of-cerebral-palsy/>)

Evidence/Guidelines/Clinical Expertise



Benefits of Supported Standing

- Hip stability
- Range of movement
- Reduce spasticity
- Improved alignment
- Bone mineral density
- Pressure relief
- Bladder and bowel
- Light physical exercise

(Paleg et al, 2013, 2021; Verschuren et al., 2014)



Standing Dosage recommendations

Standing programs 5 days per week positively affect:

- Bone mineral density (60 to 90 min per day).
- Hip stability (60 min per day in 30° to 60° of total bilateral hip abduction).
- Range of motion of hip, knee, and ankle (45 to 60 min per day).
- Spasticity (30 to 45 min per day).

- Very specific goals focussed only on the body.
- No maximum amount of time identified in literature
- How often do we prescribe '60 minutes day, 5 times per week'?
- Authors also advise promotion of Activity, Participation, Social Inclusion

(Paleg et al, 2013)

Expert Consensus/International Guidelines for Cerebral Palsy – Hip Health

22	Musculoskeletal	Conditional	For	Regular use of standing equipment for positioning	We suggest regular use of standing equipment for positioning as part of an active intervention program, to potentially decrease hip migration percentage and maintain hip abduction range of motion. ⁸⁹⁻⁹¹	Children with or at high risk of CP and risk of hip dysplasia, at the age when weight-bearing in standing would begin in children who are non-weight bearing
----	-----------------	-------------	-----	---	---	--

(Morgan et al, 2021)

Age in months	1-3	3-6	9	12	18	24	36	48	Across the Lifespan
Postural Management in Lying									
Adaptive Seating									
Standing Device									
Supported Stepping									
Bathing/Toileting									
Lift Systems									

Fig. 1. Recommended schedule of interventions.

(Paleg & Livingstone, 2022)

Holistic benefits of supported standing

- Development of skills
- Communication, respiratory, visual (Goodwin et al., 2018)
- Facilitate activity and participation (McClean et al, 2022)
- Exploration of environment
- Play
- Sensory feedback
- Learning
- Interaction and socialisation
- Fun/enjoyment



Which Standing Frame do we choose?



Prone/Upright Standers

- Can help develop head and trunk control
- A more independent transfer into equipment
- Lower limb and trunk stretching
- More optimal amount of weightbearing through lower limbs
- Prone angle can be useful for developing upper trunk and head control
- Participation



Supine Standers

- Usually provides more head and upper trunk support
- Requires a mechanical hoist for transfer as child gets older
- Often more features for deformity/contracture support
- A slower progression from into standing if required
- Generally a more supportive way of standing: GMFCS level V



Multi-position Standers

- Allow for multiple positional changes
- Alternative exploration/viewpoint of environment
- Allow for varying ways to participate
- Alternative proprioceptive feedback
- Options for progression in ability (or regression)
- Potentially support multi-faceted goals better
- Could be suitable for multi-use environments



Individualised Prescription



Copyright Etac, retrieved <https://www.etac.com/en-gb/uk/products/paediatrics/standing/r82-meerkat/>



Copyright Made for Movement Group.
Retrieved from Schmidt-Luke et al, 2018



Copyright Baffin Retrieved from
https://www.shop.baffin.co.uk/seating-systems/23-baffin-trio-.html?gclid=EAlaIqobChMln8UJ8rG_QIVA8btCh1o1wg-EAQYASABEgLIy_D_BwE



Copy right Berollka, retrieved from
<https://berollka.de/en/produkt/physio-big-relax-motion/>




Copyright Ornessa, retrieved from
https://www.ormesa.com/en/_app/



Starting the Journey



The F Words



Indications for starting supported standing in infants

- Evident gross motor delay and predicted delay in standing
- Unable to stand independently at typical developmental window
- Limited weightbearing through lower limbs
- Hip stability considerations
- Tone management and range of movement maintenance
- 24-hour postural care plan
- Skill development
- To encourage participation and inclusion with peers and family



Precautions and Contraindications

- Consent not given
- Fractures (unless medical clearance has been given)
- Medically unstable
- Significant pain
- Extreme fatigue
- Behaviour and Safety
- Postural hypotension
- Seizure Activity
- Contractures
- Hip instability
- Post Orthopaedic Surgery (Medical clearance and post op instructions)

Factors to consider

- Age and gross motor milestones
- Head control
- Trunk Control
- Current transfer ability
- Orthotics
- Potential development and prognosis
- Behaviour
- Carer and child goals
- Utilising guidelines holistically
- Growth
- Where is it going to be used?
- Who is going to use it?
- Other equipment within the carer environment
- F-Words!




MDT Collaboration



The Significance of Communication

- Have we jointly set goals with the family/carers?
- Are the goals child-centered?
- Any focus on participation, activity and home/nursery?
- Have we demonstrated the benefits of the standing system that AREN'T focused on body structure/impairment?
- Carer's understanding of why you are trying a standing frame
- Listen!!
- Empathise!




Begin Gradually


- Take time to build tolerance, understanding, routine, engagement, assess impact on body.
- Carers need time, education and support.
- Assessment and prescription is not always perfect first time! Don't be scared to re-evaluate and re-assess.
- Ultimate dosage should consider Body Structure and Function/Fitness goals, but most important goals to child and carer are most often Activity and Participation related.

Make it fun!

- Do we routinely just assess and prescribe the standing frame for static posture?
- Do we use a standing frame as part of therapy sessions?
- Do we show how it can be used for play, exploration, interaction?



Carers/Family have the best ideas!



Jenx and Jiraffe Standing Assessment Toolkit



Standing Frame Assessment Toolkit



Standing Frame Assessment Form

Any recent investigations?:

Medication and treatment:

Child/carer expectations or concerns:

Key Questions

	Yes	No
Any current fractures?	☐	☐
Any history of fractures?	☐	☐
Known osteoporosis?	☐	☐
Seizures?	☐	☐
History of extreme fatigue?	☐	☐
History of pain?	☐	☐
Any history of fainting?	☐	☐
Postural blood pressure concerns?	☐	☐

	Yes	No
Lower limb or spinal orthotics?	☐	☐
Severe lower limb contractures?	☐	☐
Pressure care concerns or injuries?	☐	☐
Currently unwell?	☐	☐
Hip stability concerns?	☐	☐
Any recent surgery? (e.g. Adductor Pexesis, Tendonitis, DOL, Gastrostomy, Tracheostomy)	☐	☐

Standing Frame Assessment Form

Sensation and sensory needs:

Behaviour:

Current transfer ability:

Environment:

Any previous standing frame use? supine/prone/sit/position/dynamic etc

If yes, known type: i.e. brand,

Yes ☐ No ☐

Successes/difficulties with previous standing frame, and what led to re-assessment?:



References

- Beopoulos A et al. Autism spectrum disorders pathogenesis: Toward a comprehensive model based on neuroanatomic and neurodevelopment considerations. *Front Neurosci.* 2022; 3;16:988735. <https://pubmed.ncbi.nlm.nih.gov/36408388/>
- Bradbury, M. et al. "Tools for participation: living aids and the F-words for childhood development." *Paediatrics and Child Health.* 2021; 31(9): 352-358. <https://www.sciencedirect.com/science/article/abs/pii/S1751722221000962>
- Bull F.C. et al. World Health Organization 2020 guidelines on physical activity and sedentary behaviour. *British Journal of Sports Medicine.* 2020; 54: 1451-1462. <https://pubmed.ncbi.nlm.nih.gov/33239350/>
- Carbone, P. S. et al. Promoting the Participation of Children and Adolescents With Disabilities in Sports, Recreation, and Physical Activity. *Pediatrics.* 2021; 148 (6): e2021054664. <https://pubmed.ncbi.nlm.nih.gov/34851421/>
- Downs, J., et al. Environmental enrichment intervention for Rett syndrome: an individually randomised stepped wedge trial. *Orphanet J Rare Dis.* 2018; 13, 3: 1-9. <https://ojrd.biomedcentral.com/articles/10.1186/s13023-017-0752-8#citeas>
- Gericke, T. Postural Management for children with cerebral palsy: consensus statement, *Dev Med Child Neurol.* 2006; 48(4): 244. <https://pubmed.ncbi.nlm.nih.gov/16542509/>
- Goodwin, J, et al. Understanding frames: A UK survey of parents and professionals regarding the use of standing frames for children with cerebral palsy. *Child: care, health and development.* 2018, 44(2), 195–202. <https://doi.org/10.1111/cch.12505>
- Goodwin J., et al. Understanding frames: A qualitative exploration of standing frame use for young people with cerebral palsy in educational settings. *Child Care Health Dev.* 2019; 45(3):433-439 <https://pubmed.ncbi.nlm.nih.gov/30883842/>
- Harvard University Center on the Developing Child. *Brain Architecture.* 2023. retrieved from <https://developingchild.harvard.edu/science/key-concepts/brain-architecture/>.
- McLean, LJ, Paleg, GS, Livingstone, RW. Supported-standing interventions for children and young adults with non-ambulant cerebral palsy: A scoping review. *Dev Med Child Neurol.* 2022; 00: 1– 19. <https://onlinelibrary.wiley.com/doi/epdf/10.1111/dmcn.15435>

References

- Morgan, C., et al. Effectiveness of motor interventions in infants with cerebral palsy: a systematic review. Dev Med Child Neurol, 2016; 58: 900-909. <https://onlinelibrary.wiley.com/doi/10.1111/dmcn.13105>
- Morgan, C. et al. Early Intervention for Children Aged 0 to 2 Years With or at High Risk of Cerebral Palsy: International Clinical Practice Guideline Based on Systematic Reviews. JAMA Paediatrics. 2021; 175(8), 846-858. <https://pubmed.ncbi.nlm.nih.gov/33999106/>
- NICE. Spasticity in Children and Young People with Non-Progressive Brain Disorders: Management of Spasticity and Co-Existing Motor Disorders and Their Early Musculoskeletal Complications. 2012. (NICE Clinical Guidelines, No. 145.). London: RCOG press. <https://www.nice.org.uk/guidance/cg145>
- NICE. Cerebral Palsy in under 25's: Assessment and Management. 2016. <https://www.nice.org.uk/guidance/ng62>
- Novak I. Evidence-based diagnosis, health care, and rehabilitation for children with cerebral palsy. J Child Neurol. 2014 Aug;29(8):1141-56. <https://pubmed.ncbi.nlm.nih.gov/24958005/>
- Novak, I. et al. State of the Evidence Traffic Lights 2019: Systematic review of interventions for preventing and treating children with cerebral palsy. Current Neurology and Neuroscience Reports. 2020; 20:3. <https://pubmed.ncbi.nlm.nih.gov/32086598/>
- Novak I. et al. Early, Accurate Diagnosis and Early Intervention in Cerebral Palsy: Advances in Diagnosis and Treatment. JAMA Pediatr. 2017;171(9):897-907. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC9641643/>

References

- Paleg, G., Smith, B. & Glickman, L.. Systematic review and evidence-based clinical recommendations for dosing of pediatric supported standing programs. *Pediatr Phys Ther.* 2013; 25(3):232-47.
https://www.researchgate.net/publication/241697453_Systematic_Review_and_Evidence-Based_Clinical_Recommendations_for_Dosing_of_Pediatric_Supported_Standing_Programs
- Paleg, G. & Livingstone, R. Evidence-informed clinical perspectives on postural management for hip health in children and adults with non-ambulant cerebral palsy. *Journal of Pediatric Rehabilitation Medicine: An Interdisciplinary Approach Throughout the Lifespan.* 2022; 15(1): 39–48. <https://pubmed.ncbi.nlm.nih.gov/35275575/>
- Park, H.P. et al. Sedentary Lifestyle: Overview of Updated Evidence of Potential Health Risks. *Korean J Fam Med.* 2020; 41(6): 365-373. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7700832/>
- Smith, B. et al. Physical activity for general health in disabled children and disabled young people: summary of a rapid evidence review for the UK Chief Medical Officers' update of the physical activity guidelines. Department of Health and Social Care: London, UK. 2022. Available at <https://www.gov.uk/government/publications/physical-activity-in-disabled-children-and-disabled-young-people-evidence-review/physical-activity-for-general-health-benefits-in-disabled-children-and-disabled-young-people-rapid-evidence-review>
- Vance. The Emerging Role of Cognitive Remediation Therapy, Activities, Adaptation & Aging. 2009; 33:1, 17-30.
<https://www.tandfonline.com/doi/abs/10.1080/01924780902718541>
- Verschuren, O. et al. Muscle activation and energy-requirements for varying postures in children and adolescents with cerebral palsy. *J Pediatr.* 2014; Nov; 165(5):1011-6. <https://pubmed.ncbi.nlm.nih.gov/25151195/>
- World Health Organization. 2019. WHO guidelines on physical activity, sedentary behaviour and sleep for children under 5 years of age. Available at: <https://www.who.int/publications-detail/guidelines-on-physical-activity-sedentary-behaviour-and-sleep-for-children-under-5-years-of-age>

Thank you